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Type E 1500, E 2000, E 2800, E 4000, E 6000, E 7000, E 11000, E 15000

Read these operating instructions carefully before installing and starting up the system!

Kadotec Wassertechnik part of AQUA Management GmbH
Emil-von-Behring-Str. 3, 63128 Dietzenbach
Tel. +49 (0) 6074 8053113 * Fax +49 (0) 6074 8053111 *
E-Mail info@kadotec.de
www.kadotec.de

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Preface

Dear Sir or Madam,

The Kadotec team thanks you for the confidence you have placed in us.

With your selection of this mixed bed ion exchanger system you have chosen an innovative and high-quality product that will provide a long service life.

Please study the information that is given in these operating instructions on installation and operation before you start to install and start this water purification system.

Please also only appoint informed and trained personnel to operate the system and only use Kadotec original accessories, replacement parts and consumables.

The above requests are particularly important as we, the manufacturer, cannot accept liability for any damages to the instrument, connected peripherals, building or persons that occur because of improper installation, operation or use of the system, or of the use of parts of any type of other manufacturers.

We wish you every success with and because of your new water purification system!

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
1 Contents

2	Parts list / options	5	18	Installation and operation	12
13	Symbols used	5	18.1	Putting into operation	12
14	Safety precautions	10	19	Regeneration	12
15	Intended use	11	20	Specifications E1500 – E 6000	13 21
16	Process description	11		Specifications E 7000 – E15000	13 22
17	Determination of the exchange capacity	12		Accessories for ion exchangers	13 24
				EU – Declaration of Conformity	14 25
				Notes	15

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Parts list

- ① Raw water connector
- ② Raw water hose R ¾"
- ③ Pure water hose R ¾"
- ④ Measuring cell with R ¾" T-connector
- ⑤ Kadotec DC 400 conductivity meter
- ⑥ Transformer plug
(only in scope of delivery art-no. 14180500)
- ⑦ Raw water solenoid valve with R ¾" connection
(only in scope of delivery art-no. 14180700)
- ⑧ Power plug with switch socket for solenoid valve
connection
(only in scope of delivery art-no. 14180700)
- ⑨ Power plug with potential free contact,
RS 232 + USB
(only in scope of delivery art-no. 14160600)
- ⑩ Kadotec analogue conductivity meter
(only in scope of delivery art-no. 14160100)
- 1 Gasket R ¾"
(only in scope of delivery art-no. 14160100)
- 1 Mixed-bed ion exchanger
E 1500 – E 15000

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13 Symbols used

This warning symbol is given for a danger that, when not avoided, could result in slight or moderate injury.

This warning symbol is given for a danger that, when not avoided, could result in slight or moderate injury.

This symbol refers to particularly useful notes.

14 Safety precautions

- Heed these notes and precautions for your own well-being and proper functioning of the cartridge!
- Read these operating instructions carefully before installing and starting up the system and pay particular attention to the safety precautions.
- Improper use of this system or use of it in any other way than that described in this user manual invalidates the guarantee.
- When a constructional change is made to the system, or a product of a different manufacture is installed in it, the CE-mark is invalidated.
- The intended purpose of mixed-bed ion exchangers E 1500 – E1 5000 is to deionize potable water.
- 30° C is the maximum permissible temperature for the potable water used as feedwater for mixed-bed ion exchangers DS 1500 – DS 15000. The guarantee does not cover damage to joints, gaskets and/or ion-exchanger material caused by water of a higher temperature (leakages).
-  The room in which the water deionizer is installed must be frost-free and the temperature in it is never to go below + 2° C.
-  The installation area must have a floor drain with at least DN 50 pipe. When no floor drain is available, a water watcher (article no.: 13012900) is to be installed. The manufacturer is not liable for damages caused by overflow of water.
- For safety reasons, shut the water supply off when the system is not required for a longer period (weekends or holidays).
-  Do not open the central cover lid (D) in pressure conditions on the cartridge. The opening should only be carried out by an authorized regeneration dealer.
-  For safety reasons it must be a back flow preventer and pipe cut-off combination to prevent return flow to the public potable water supply, as is required by the DIN EN 1717 part 4 regulations for potable water installations.
- Follow all currently appropriate laws and regulations, including those on accident prevention, which are in force at the location where the system is to be installed.
-  The installation of a particle filter in the ion exchanger outlet is recommended to protect downstream users from possible harm.
-  Under some circumstances, a defective cartridge could allow ion exchange resin to leak out.
CAUTION: Danger of slipping on leaked-out resin!
As the resin material can cause slight irritation on contact with skin and/or eyes, immediately turn off the raw water tap to stop further leakage and proceed further as follows:
On skin contact, immediately wash contacted skin areas with plenty of water.
On contact with eyes, remove contact lenses if this is possible, then immediately and carefully wash eyes with plenty of water for about 10 minutes. When irritation still occurs, seek medical advice immediately.
Wear protective gloves to collect up the leaked out resin, fill it into a clearly labelled plastic container and return it to the authorized specialist supplier for waste disposal.
- The deionized water that is produced is not fit for drinking.

15 Intended use

The mixed bed ion exchanger is a water treatment system for the purification of drinking water to deionized water, especially for use for:

- Dispensing chemist's
- Printeries
- Laboratories
- Laboratory washing machines
- Humidifiers
- Medical technology
- Hospitals
- Battery water producers
- Reverse osmosis product water
- Pharmaceutical industry
- Industrial cleaning plants
- Steam generators
- etc.

The water that is produced is not fit for human consumption.

This system is not a medical device.

16 Process description

Kadotec mixed bed ion exchangers work according to the ion exchange principle. The dissolved salts that are present in the water are in the form of ions. The ion exchange resins exchange these for H^+ and OH^- ions. The result is pure, demineralized water. The resins used have a certain ion exchange capacity, whereby the volume of raw water that can be processed depends on the total salt content of the

raw water. A conductivity meter that is fitted to the outlet of the cartridge containing the resins shows the degree of resin exhaustion. When the reading reaches the limiting value, the ion exchange capacity has been reached. The cartridge containing the ion exchange resins is then to be sent to a central regeneration station for appropriate regeneration treatment and return.

17 Determination of the exchange capacity

Feed water quality

With a feed-water quality of $1^\circ dH = 1.25 eH$ total dissolved solids (approx. $30 \mu S/cm$), the exchange capacity is ten times the type number of the particular ion exchanger type.

Using E 4000 as example:

- E4000 at $30 \mu S/cm$ (approx. $1.25^\circ eH$ total dissolved solids) = 40 000 liters
- E 4000 at $300 \mu S/cm$ (approx. $12.5^\circ eH$ total dissolved solids) = 4000 liters
- E 4000 at $600 \mu S/cm$ (approx. $25^\circ eH$ total dissolved solids) = 2000 liters

Please note:

The exchange capacity is lowered by high amounts of carbonic acid and silicic acid in the feedwater.

Pure water quality

The pure water taken from the ion exchanger is free of salts and carbon dioxide. It therefore has a pH of approximately pH 7. When a lower pH is measured in the pure water, this is usually the result of measurement error.

The reason for this is that pure water has a craving for carbon dioxide and can take up large amounts of it. Even when it only takes up extremely low amounts, the pH is lowered. Only flow-through pH measurement procedures that are carried out under the exclusion of air bring correct results. It is clear from this that pH measurements with indicator paper do not provide reliable results.

18 Installation and operation

Site prerequisites:

Install the mixed bed ion exchange cartridge so that it is rigidly held. Ensure that it is freely accessible.

1. The installation area must be frost-free, minimum temperature + 2° C
2. Potable water stopcock with R ¾" male thread
3. DN 50 waste pipe or water watcher with sensor
4. No source of heat in the vicinity of the ion exchange system or peripheral equipment
5. Use only the hoses supplied. The use of materials that are not suitable for pure water, (e.g. nonferrous metal) could possibly lead to deterioration over time.
6. The cartridges are solely for connection to a potable water supply.

Connection:

7. Remove the cartridge from the cardboard box and save the box for return of the exhausted cartridge for regeneration.
8. Remove the closing caps from the cartridge connectors and save them for fitting them to the exhausted cartridge.
9. Fit the dosing orifice (only supplied in the additional operating instructions consignment for cartridges E 1500 to E 4000) in the hose connecting nipple of the raw water hose.

CAUTION: Proper system functioning and the high quality of the pure water can only be attained with inserted dosing orifice position.

10. Screw the connecting nipple of the raw-water hose with inserted dosing orifice to the potable water shut off valve.
11. Connect the free end of the raw-water hose to the connector labelled **INLET (A)** of the mixed-bed ion exchange cartridge. Only hand tighten the hose nipple.
- 11a. When the conductivity meter has a built-in measuring cell: Screw the conductivity meter with built-in measuring cell to the cartridge connector labelled **OUTLET (B)**. Connect the free connector of the conductivity meter with the pure water hose that is supplied.

- 11b. When the conductivity meter has a separate measuring cell: Screw the measuring cell to the connector of the mixed-bed ion exchanger labelled **OUTLET (B)**. Connect the free connector of the measuring cell to the pure water hose that is supplied. Only hand-tighten the hose nipple.
12. Connect the free end of the pure water hose to your user.

CAUTION: Do not mix up the raw water hose **INLET (A)** with the pure water hose **OUTLET (B)** as the high quality of the product water would then not be attainable.

18.1 Putting into operation

Air-vent the cartridge:

13. Open the **vent-valve (C)** of the central **cover lid (D)** of the cartridge
14. Carefully open the raw water tap and wait until water emerges from the vent-valve.

CAUTION: Do not screw the **vent valve (C)** out of the thread but just turn it anti-clockwise until air emerges.

15. When water emergencies, close the vent-valve.
16. The mixed-bed ion exchanger is now ready for use.
17. The quality of the water that is produced can be read off from the conductivity meter.

CAUTION: The values of the conductivity are irrelevant when no product water is drawn out. True values are only shown when water flows through the measuring cell.

18. The mixed-bed ion exchanger cartridge must be regenerated when the maximum conductivity value for the particular application is reached.

19 Regeneration

We recommend to keep an additional cartridge in reserve to continually have pure water at your availability.

To maintain the quality of the ion exchangers, the regeneration is carried out by our authorized dealers with stakpure mixed bed resin. Upon request, we will gladly inform you of the dealer responsible for you.

20 Specifications E 1500 – E 6000

Kadotec mixed bed ion exchanger	E 1500	E 2000	E 2800 pressure resistant	E 4000	E 6000
Water demineralizer			pressure resistant up to 10 bar		
Stainless-steel cartridge					
Flow rate	300 l/h	300 l/h	950 l/h	1000 l/h	1000 l/h
Capacity at 10°dH total dissolved solids	1500 l	2000 l	2800 l	4000 l	6000 l
Pure water quality			0.1 – 20 µS/cm		
Max. water temperature			30° C		
Max. Inlet pressure			6 bar		
Material			Stainless steel – AISI 316		
Dimensions mm (Ø x height)	240 x 410	240 x 490	240 x 570	240 x 700	240 x 1155
Connection*			R ¾"		
Weight, approx.	14 kg	18 kg	24 kg	27 kg	45 kg
Order-No.	12115000	12200000	12280000	12400000	12600000
Dosing orifice	31000013	31000013	31000032	31000033	n/a
Wall mount			13140500		

*On request, the cartridges can also be equipped with quick-connect couplings.

21 Specifications E 7000 – E 15000

Kadotec mixed bed ion exchanger	E 7000	E 11000	E 15000
Water demineralizer		pressure resistant	
Stainless-steel cartridge		pressure resistant up to 10 bar	
Flow rate	2000 l/h	2500 l/h	3000 l/h
Capacity at 10°dH total dissolved solids	7000 l	11000 l	15000 l
Pure water quality		0.1 – 20 µS/cm	
Max. water temperature		30° C	
Max. Inlet pressure		6 bar	
Material		Stainless steel – AISI 316	
Dimensions mm (Ø x height)	363 x 600	363 x 850	363 x 1100
Connection*		R ¾"	
Weight, approx.	55 kg	70 kg	90 kg
Order-No.	12700000	12110000	12150000

*On request, the cartridges can also be equipped with quick-connect couplings.

22 Accessories for ion exchangers

Designation	Order no.
Conversion set to quick connect nipple, nickel-plated	13150000
Conversion set to quick connect coupling hose, nickel-plated	13150600
Conversion set to R ¾" quick connect nipple with check valve	13150700
Conversion set to quick connect coupling – hose with check valve	13150800
2 x 1.5m hose, DN 10 R ¾"	13169000
2 x 1.5m reinforced hose, DN 20 R ¾"	13169100
Wall bracket for cartridge types E 1500 –E 4000	13140500
Water watcher with solenoid valve and sensor R ¾"	13012900
Pressure reducer, 0 – 4 bar R ¾"	13141400
Solenoid valve, nickel plated R ¾"	13150300
Pure water tap, straight nozzle R ¾"	13140000
Pure water tap, angled nozzle R ¾"	13140100
Pure water shut off valve R ¾"	13140300
Pure water two-way distributor R ¾"	13142000
Pure water three-way distributor R ¾"	13140200

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24 EU-Declaration of Conformity

in the sense of **directive 2014/68/EU** for pressure vessels.

Manufacturer: Kadotec Wassertechnik Part of AQUA Management GmbH, Emil-Von-Behring-Str. 3, 63128 Dietzenbach

We hereby certify that the conception and form of the machine named below that we have brought to market are in accordance the determinations of the named article of the EC directive. This declaration is no longer valid when changes are made to the product.

Product designation: Mixed bed ion exchanger
Type: E 1500, E 2000, E 2000-17, E 2800, E 2800-22, E 4000, E 6000, E 7000,
E 11000, E15000
Art. No.: 12115000, 12200000, 12200102, 12280002, 12280102, 12400000, 12600000,
12700000, 12110000, 12150000

The following harmonised standards were applied:
DIN EN 13445-3:2018-12

The following national standards, directives and specifications were applied:
AD 2000 set of regulations

Dietzenbach, 10.01.2022

Karlheinz Freihoff, Executive Director

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Hersteller: Kadotec Wassertechnik Part of AQUA Management GmbH, Emil-Von-Behring-Str. 3, 63128 Dietzenbach